



**S**  
**SIALOR**  
M I L A N O

**YALODÉ FILL**  
Sialor Milano for aesthetic medicine



**SIALOR S.r.l.**

Via Lodovico Mancini, 5  
Milano - ITALY  
info@sialormilano.it  
www.sialormilano.it




**IT BECOMES  
ESSENTIAL TO REDEFINE  
THE CONTOUR OF  
THE FACE WHILE  
RESPECTING  
ITS NATURAL HARMONY**

In young skin, the dense network of Collagen, continually renews itself giving support and structure to the tissues.

**Around age 40, however, the oval of the face begins to be affected by gravity.**

The skin sags, the cheeks become hollow, and the lower part of the face thickens. The result is a tired face with a less defined contour and wrinkles and furrows of expression.



## INTERPRETING PATIENTS' DESIRES IS THE MISSION OF SIALOR MILANO

The search for Harmony to enhance one's Beauty requires personalized and exclusive treatments.

The woman desires well-delineated facial contouring, radiant skin and youthful, fresh appearance.

The man may be interested in well-delineated masculine profile and well-defined jawline.

The woman and man want to appear natural and never artificial.


## YALODÉ FILL

### EXCLUSIVE FILLER WITH TRIPLE EFFECT

Yalodé Fill offers one of the most High Molecular Weight 2,000 kDa biotechnology-derived Linear Hyaluronic Acid formulations on the market.

Yalodé Fill treats physiological depletion of **Hyaluronic Acid at the dermal level and with its High Molecular Weight** works to improve skin texture and laxity.

It adds gentle volume by redefining contours that have lost their natural harmony.



## THE EFFECTIVENESS OF YALODÉ FILL

The hygroscopic and viscoelastic properties of High Molecular Weight Hyaluronic Acid allow the hydration of tissues favoring their remodeling with corrective effect of tonicity and scarring outcomes.


It is a very comprehensive anti-aging system that ensures the turgor of the dermis through the progressive improvement of skin elasticity.

**YALODÉ FILL is intended for the treatment of skin imperfections and particularly in specific areas of the face.**

- Nasolabial folds
- Peribuccal wrinkles
- Glabellar wrinkles
- Frontal wrinkles
- Facial scars including acne scars
- Areas of the face in need of facial tissue enrichment (cheeks, chin, cheekbones, lips) with temporary volumetric augmentation by soft tissue augmentation.




## THE ACTION OF YALODÉ FILL



EYE CONTOUR  
WRINKLES



NASOLABIAL  
FOLDS



NECK AND FACIAL  
CONTOUR  
SKIN LAXITY



## INJECTIVE TECHNIQUE DEDICATED TO EACH FACE

**FOR A PATIENT TO BE SATISFIED, IT IS  
NECESSARY TO UNDERSTAND THEIR  
EXPECTATIONS.**

With Yalodé Fill it is possible to personalize the treatment while respecting the naturalness of each face.



### INJECTION DEPTH

The implant should be placed in the Mid-Surface Derma using a 30 G x 13mm needle.

### GOOD MOLDABILITY

The cohesive gel matrix is easily adaptable and moldable. Light massage can be implemented to help distribute even better the injected product.

**VISIBLE RESULTS**  
with an extended duration  
of up to 12 months

**TRIPLE EFFECT**  
vertical lifting, contouring,  
neocollagenesis

**HIGH ELASTICITY AND VISCOSITY**  
lifting and contouring action



YALODÉ FILL 2% - YALODÉ FILL 2.5%



**Prevents the reduction of Hyaluronic acid**, promoting BIORISTRUCTURE of skin tissue.

**Safe and well tolerated**, improves hydration, firmness and skin elasticity.

**Available in 2 different concentrations** of cross-linked Hyaluronic acid to ensure treatment in every therapeutic need.

**ERGONOMIC DESIGN**  
for enhanced handling  
by the clinician.


**CONTROL OF  
INJECTED MATERIAL**  
by graduation on the  
syringe.

**AGO STABILITY**  
by the perfect seal  
of the Luer -Lock  
attachment.



YALODÈ FILL **2%**      YALODÈ FILL **2.5%**

<b>Indications for use</b>	Lip contour wrinkle treatment	Treatment of facial wrinkles
<b>Concentration</b>	20 mg/ml	25 mg/ml
<b>Injection zone</b>	Derma preferably the deeper layer	Deep supraperiosteal hypodermis
<b>pH</b>	6,00 – 7,50	6,00 – 7,50
<b>Needle</b>	27G 1/2"	27G 1/2"
<b>Package</b>	1 x 1 ml	1 x 1 ml
<b>Duration of Result</b>	4 - 8 months	6 - 10 months
<b>Conservation</b>	5 – 25° C	5 – 25° C





## THE CORRECTIVE EFFECT OF YALODÈ FILL (2.-2.5%) IS IMMEDIATELY VISIBLE AFTER TREATMENT

The rehabilitative action effectively corrects the emptying cheeks syndrome in the stands of naso-labial wrinkles, Characterizes the chin and harmonizes the oval of the face, redefining the jaw line.

Yalodè Fill (2-2.5%) acts mechanically by absorbing water generating swelling and biologically by increasing the synthesis of neoCollagen and Glycosaminoglycans.

The dual effect makes it possible to generate volume that immediately and permanently corrects skin depressions.



## SAFE & EFFECTIVE

### HIGH QUALITY STANDARDS

- **Cross-linked Hyaluronic Acid** of biotechnological origin, definitely not animal.
- **Maximally purified cross-linked Hyaluronic Acid.**

<b>Nucleic acids</b>	≤ 0,2%
<b>Proteins</b>	≤ 0,1%
<b>Chlorides</b>	≤ 0,03%
<b>Iron</b>	≤ 80 ppm
<b>Microbial contamination (TAMC)</b>	≤ 100 CFU/g
<b>Bacterial endotoxins</b>	≤ 0,05 I.U./mg

### ASSURED STERILITY OF PACKAGING

- Low risk of allergic reactions.
- Sterilization by steam and irradiation.
- Increased durability of the implant.



## SUPPORTING TECHNOLOGY

The half-life of cutaneous Hyaluronic acid does not exceed 48 hours in equilibrium between production and demolition by hyaluronidases.

To overcome this limitation, in aesthetic medicine, Hyaluronic acid, for intradermal use, is modified by cross-linking **to obtain a product that is more resistant to demolition and therefore more durable: cross-linking.**

The activity of cross-linked Hyaluronic acid, not only due to the function of restoring volume loss and filling wrinkles, but also to a protracted stimulation that, after time, allows further improvement of the treated skin tissue and therefore constitutes a plus of considerable importance.

## CROSS-LINKING TECHNOLOGY

Cross-linked Hyaluronic acid has bridges connecting HA molecules, with the formation of a three-dimensional structure. This system exhibits a strong resistance against the degradative processes of Hyaluronidases, and consequently greater durability.

BDDE is used as the cross-linking agent, which forms strong covalent bonds between the hydroxyl groups of the HA chains, increasing the strength and rigidity of the 3D network.

## UNIFORM GRAIN SIZE

HA is scrupulously selected at the granulometric level to ensure a well-defined and consistent average particle size to assure a product **characterized by a high quality standard.**

# MBP

## MULTI-ACTIVE BEAUTY PATH IMPROVE AND PROLONG THE RESULTS OF AESTHETIC MEDICINE

### YALODÉ INTENSIVE WRINKLE REPAIR

Consistent use of **YALODÉ INTENSIVE WRINKLE REPAIR**, with 19 Active Ingredients, moisturizing, restorative and long-lasting with a botulinum-like effect, can promote the smoothing of wrinkles and expression lines, visibly improving the firmness, radiance and elasticity of the skin of the face, neck and décolleté.

**HYALURONIC ACID** with 3 different molecular weights, obtained through biofermentative technology, acts at the different layers of the skin.

**7 BIOMIMETIC PEPTIDS** chosen from the most studied and tested active ingredients to effectively fight wrinkles.

**CERAMIDES** promote skin hydration.

**CENTER TRITERPENES** promote Collagen synthesis and reduce the depth of wrinkles.

**MARINE COLLAGEN** provides functional amino acids for neoCollagen formation.

**DIMETHYLAMINOETHANOL** natural antioxidant that tones muscle making skin firmer and more elastic.

**BATYL ALCOHOL** protects the skin from damage caused by UV radiation and external aggressions.

**4 VEGETABLE SUBSTANCES** Bisabolol, Shea Butter, Jojoba Oil and Grapeseed Oil with specific activities.

### TRIPLE LAYERS TABLETS

**YALODÉ TRIPLE LAYERS TABLETS** provides from within, where creams don't reach, targeted Natural Nutrients that combine the super-hydrating and anti-aging action of Hyaluronic Acid with cutting-edge Functional Actives and Antioxidants to promote skin firmness, radiance and elasticity to protect cells from oxidative stress.

**Serving size: 1 tablet daily with water**

High molecular weight **HYALURONIC ACID** of biofermentative origin.

**RED GRAPE VINE** plant species with antioxidant function whose seeds contain Proanthocyanidins.

**RESVERATROL** natural polyphenol, is the major component of Polygonum Cuspidatum, a plant substance with antioxidant action.

**COENZYME Q10** naturally present in all cells.

**ORGANIC SILICON** constitutes the readily available source of Silicon for humans.

**BIOTIN** and **ZINC** help maintain normal skin tone.

**VITAMIN C** useful for the formation of Collagen.

**SELENIUM** and **VITAMIN E** protect cells from oxidative stress.



**Hyaluronic Acid**

**Active Ingredients**

**Hyaluronic Acid**



**SIALOR S.r.l.**

Via Lodovico Mancini, 5

Milano - ITALY

[info@sialormilano.it](mailto:info@sialormilano.it)

[www.sialormilano.it](http://www.sialormilano.it)



[sialor\\_beauty](#)



[Sialor Beauty](#)



[Sialor Milano](#)